Purpose of the instructions

This Manual is relative to a elbow orthosis, in compliance with Regulation 2017/745 concerning class 1 medical devices, and intends to provide a guide for application, regulation (by orthopaedic specialist) and for safe and correct use (by patient) throughout all his life. Therefore the Declaration of Conformity has been issued by Tecnoway under its own responsibility in accordance with Annex IV of this Regulation

Identification of the manufacturer and the product

Tecnoway produces medical devices, in particular orthopaedic braces. The Firm is constantly searching for new materials innovative technologies and ergonomic solutions. Its system of quality management is certificated in accordance with ISO 9001 & ISO 13485 regulation. The product specified in this instruction for use is a functional brace for the control of elbow flexion-extension and forearm pronation-supination.

General description and mode of action

The device described in this manual is composed primarily of a frame in aluminum alloy and non-toxic material padded supports making the device light, comfortable and safe. Easily removable controls the movement of the elbow without limiting the operation of daily hygiene

Caution for use

The device should be used exclusively for orthopaedic treatment of the upper limb. The duration of daily use and the application period must be indicated in the prescription of an orthopaedic doctor. The pressure exerted by the device should not act upon injury, or caused localized symptoms caused by compression, or compress nervs fibers and blood vessels. The device is not fireproof, do not use near flames. Not use near strong electro-magnetic field. Be careful to open the Velcro fastening slowly to avoid a sudden collapse of the upper limb.

Features

Ambidextrous joint of the latest generation with patented technologies. Anticorrosive alloy bars with adjustable slide. Malleable aluminium cuffs enable the fitter to customise the Double-R.O.M. to varying anatomy. Soft padding lined with breathable fabric. Joint adjustment of flexion-extension every 10 ° with 120° range, possibility of locking of the joint. Telescopic hand support to control the pronation-supination adjustable and removable support (only Unibor 02). The orthosis stabilizes the arm and wrist ensuring a comfortable and effective therapy.

Indications

Control of flexion-extension of the elbow and prono-supination after surgery or trauma.

Application (Orthotist/patient)

1. Set the R.O.M. at 0° Flex and 120° Ext version right or left based on the limb treated acting on the cursors "C" (Fig.A).
2. Adjust the lenght of the brace and the position of the valve acting on the red buttons "D".
3. Place the brace on the arm by matching the articulation "E" with the epicondyle.
4. Adapt the aluminium valve manually around the arm
5. Close the velcro straps starting with those next the elbow.
6. Adjust the flexion-estension of the R.O.M. "E" acting on the cursors "C".
7. Unfasten the wrist strap and insert the terminal buckle of the armsling strap by velcro strap "F".
8. Slip the loop of the strap and sling around the neck, fasten the buckle "G" on the wrist.
9. If necessary adjust the lenght of the belt armsling
10. Insert the handle support in the special groove pressing the red buttons "H"(only Unibor 02)
11. Adjust the position and release the buttons
12. If necessary rotate the handle in the desired position by loosening the screw "I" (only Unibor 02)
13. Correct use is when the epicondylus supports are very adherent to the articulation and arm/forearm bands are very adherent to the muscles without create excessive pressure.

Warning

The device in question must be prescribed by an orthopaedic specialist and applied by an orthotist who will according to the needs and conditions of patient instruct on the proper and safe use of the device. Any modification or regulation different from the previously set on the device at the time if application and adjustment should be prescribed by an orthopaedic specialist and performed by an orthotist. Don't change the regulation done by the orthopaedic specialist. For hygienic and functional reason we recommend the use for single patient because it could cause risks for hygienic or functional aspect. In particular sensible person, direct contact of skin to device may cause redness and/or irritations, if those problems may occure please contact directly your doctor. If the device is no more under normal use it should be replaced and disposed off.

Incident reporting notice

Any serious incident that occurs in connection with this medical device must be reported to the manufacturer and the competent authority.

Storage instructions

When not using the device, even for short periods, store it in its original packaging to maintain a good level of cleanliness. Keep it in a cool, dry place and protect from direct sunlight.

Maintenance

Before use always check the tightening conditions of the screws, of the joints and of the adjustable elements. Check before using the state of wear of the device and its parts. If worn out contact your orthotist to substitute. Hand wash the paddings with warm water and mild soap, rinse accurately. Dry away from heat sources.

Disposal

Make sure that the device is disposed off in a correct manner.

Do not disperse in the environment.

Guarantee

The guarantee is valid only if the product is used under the conditions indicated and under the intended purposes.

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This manual has been realized by Tecnoway exclusively for its clients. Its reproduction and/or divulgation, entirely and/or partial is forbidden without Tecnoway authorization. The manual comes together with the equipment device, it must be kept in a suitable place and easily accessible. The staff assigned to application, regulation use and maintenance should follow instructions carefully all the times. Tecnoway cannot be pursued or held responsible for any damage caused by wrong application different to those prescribed in this manual.

 by **TECNOWAY SRL**
via Montefeltro n. 45 61022 Vallefoglia -PU- ITALY
+39.0721.496175 info@tecnoway.it www.tecnoway.it

COMPANY WITH SYSTEM CERTIFICATION
ISO 9001 - ISO 13485 - ISO 45001

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