**COLLARÍN ANATÓMICO CERVICAL / CONTOURED CERVICAL COLLAR
COLLIER CERVICAL ANATOMIQUE / COLAR ANATÓMICO CERVICAL
COLLARE CERVICALE BIVALVE / ANATOMICZNY KOŁNIERZ SZYJNY
ANATOMISCHE HALSWIRBELKRAUSE / ANATOMISCHE CERVICALE ORTHESE
FORMAD CERVIKAL HALSKRAGE / FORMET CERVIKAL HALSKRAGE
MUOTOILTU KAULARANGANTUKI / FORMET CERVIKAL HALSKRAVE**

961

ESP

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ENG

Fitting Instructions

FR

Instructions de pose

PT

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IT

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PL

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DE

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NL

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SV

Utprovninginstruktioner

NO

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FI

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DK

Tilpasningsinstruktion

**Selección de la talla / Size Selection / Choix de la taille
 Seleção do tamanho / Selezione della taglia
 Wybór rozmiaru / Größenauswahl / Maatkeuze
 Storleksurval / Størrelser / Kokovalikoima / Størrelser**

B

A

S 2 1/4"	25,4-33cm (10"-13")	5,46cm (2 1/4")
S 3 1/4"	25,4-33cm (10"-13")	8cm (3 1/4")
S 4 1/4"	25,4-33cm (10"-13")	10,54cm (4 1/4")
S 5 1/4"	25,4-33cm (10"-13")	13,08cm (5 1/4")
M 2 1/4"	33-40,6cm (13"-16")	5,46cm (2 1/4")
M 3 1/4"	33-40,6cm (13"-16")	8cm (3 1/4")
M 4 1/4"	33-40,6cm (13"-16")	10,54cm (4 1/4")
M 5 1/4"	33-40,6cm (13"-16")	13,08cm (5 1/4")
L 2 1/4"	40,6-48,3cm (16"-19")	5,46cm (2 1/4")
L 3 1/4"	40,6-48,3cm (16"-19")	8cm (3 1/4")
L 4 1/4"	40,6-48,3cm (16"-19")	10,54cm (4 1/4")
L 5 1/4"	40,6-48,3cm (16"-19")	13,08cm (5 1/4")
XL 2 1/4"	+48,3cm (19" Up)	5,46cm (2 1/4")
XL 3 1/4"	+48,3cm (19" Up)	8cm (3 1/4")
XL 4 1/4"	+48,3cm (19" Up)	10,54cm (4 1/4")
XL 5 1/4"	+48,3cm (19" Up)	13,08cm (5 1/4")

*Instrucciones de adaptación / Fitting instructions
Instruccions de pose / Instruções de adaptação
Istruzioni per l'applicazione / Instrukcja obsługi
Vorschrift zur npassung / Instructies voor het aanmeten
Utprovningsinstruktioner / Tilpassingsinstruksjoner Sovitusohjeet /
Tilpasningsinstuktion*

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READ CAREFULLY THESE INSTRUCTIONS BEFORE USING THE ORTHOSIS. FOR ANY DOUBT, CONSULT YOUR DOCTOR OR ORTHOPAEDIC TECHNICIAN.

Description

Contoured cervical collar with tracheotomy opening. Latex Free. Radiolucent. Easy to fit without moving the patient's neck.

Indications

In all process requiring an important neck immobilisation:

Cervical tendinitis, Osteoarthritis, Contractural diathesis, Ligamenta-muscle lesion, Medical urgencies due to cervical traumatism, Cervical post-surgery

Cautions

The doctor is the person qualified to prescribe this orthosis and decide the duration of the treatment. A healthcare professional is the qualified person to fit the product. For a correct effectiveness, you should respect the daily time of use following the doctor or orthopaedic technician indications.

For a correct function, you should use the correct size of the orthosis.

This orthosis is designed for use on one patient. Do not re-use on multiple patients.

This orthosis should never be used for any other use different to the one prescribed.

Do not allow children to play with this device.

You should inform your doctor if you suffer of some skin affections. Avoid the use of creams, in contact with the orthosis, during the treatment.

On the pressure points, the skin should not be injured nor irritated.

Check periodically that the cinch straps are tight and in good conditions (without grease, cuts, rips, etc.).

Size selection (see pag.2)

Select the size of the collar considering:

- 1.- The height measurement from the chin to the top of the sternum, in neutral position. (A).
- 2.- The neck circumference (B).

Fitting instructions (see pag.3)

This collar is composed by an anterior part (D) and a posterior one (C). Do not disassemble any other component.

Take the posterior half of the collar (C) marked "BACK ↑". Place it (C) on the patient's nape taking care not to move the patient's neck.

Then position the anterior half of the collar (D), marked "FRONT↑" and the tracheotomy opening. Place the anterior half part of the collar (D), so that the chin comes to rest in the provided recess. Always maintain the patient's neck in neutral position.

The anterior part (D) should overlap with the posterior one (C).

Adjust the width of the collar by means of the Velcro® closure so that it is adjusted with the same tension by the both sides.

Care Instructions

Clean the orthosis with a cloth, with a neutral soap and lukewarm water. Do not bleach. Dry with a dry cloth. Do not tumble dry. Do not iron. Do not dry clean.

Classification

Medical Device, Class I, non sterile, with no measurement function. This product meets the requirements of the Medical Device Directive 93/42/CEE, by which the medical products are regulated.

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Prim, S.A.
PRIM, S.A.