



27th February 2025

To Whom it May Concern

After conducting internal wash trials on Isotoner & Jura compression gloves, I can confirm that the gloves can be machine washed in a gentle cycle at 40°C.

However, we do suggest that, for environmental reasons, a cooler 30°C wash is used.

Fabric Conditioners and Bleaching Agents should be avoided.

DG53022-0	Isotoner Gloves X-Small Full Finger
DG53022-1	Isotoner Gloves Small Full Finger
DG53022-2	Isotoner Gloves Medium Full Finger
DG53022-3	Isotoner Gloves Large Full Finger
DG53023-0	Isotoner Gloves X-Small Full Finger
DG53023-1	Isotoner Gloves Small Full Finger
DG53023-2	Isotoner Gloves Medium Full Finger
DG53023-3	Isotoner Gloves Large Full Finger
82192	Thermoskin Thermal Compression Gloves X Small
83192	Thermoskin Thermal Compression Gloves Small
84192	Thermoskin Thermal Compression Gloves Medium
85192	Thermoskin Thermal Compression Gloves Large
86192	Thermoskin Thermal Compression Gloves X large
87192	Thermoskin Thermal Compression Gloves XX Large

Kind Regards,



Mr Kevin Birt
Managing Director

